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<div><span>en</span></div>
----------------------------

**Composition**
White to greyish white powder composed of tricalcium silicate, zirconium oxide and povidone.

Clear and colourless solution of calcium chloride and polycarboxylate.

**Intended use**

Permanent root canal sealer

**Clinical benefit**

Preservation of the tooth functionality

**Intended population**

Permanent and mature teeth

**Indications for use**
Permanent root canal filling in combination with gutta-percha points in case of inflamed or necrotic pulp.
Permanent root canal filling in combination with gutta-percha points following a retreatment procedure.
BioRoot™ RCS is suitable for use in single cone technique or cold lateral condensation.

**Contra indication**
Hypersensitivity to one of components
Immature or temporary teeth

**Performance characteristics**

**Technical performance**

Seal hermetically to the root canal.

**Technical specifications**

Setting time from 45 minutes to 4 hours allowing a sufficient working time.

Radiopaque, non-staining, dimensionally stable, non-soluble in fluids, adhesion to mineralized tissues and gutta-percha.

**Propertes**

- BioRoot™ RCS is highly biocompatible and reduces the risk of adverse tissue reaction.
- BioRoot™ RCS crystallization creates a tight seal within the dentin tubules.
- BioRoot™ RCS fosters mineralization of the dentin structure and facilitates the physiologic bone process. BioRoot™ RCS therefore creates a favorable environment for periapical healing.
- BioRoot™ RCS limits bacterial propagation.
- BioRoot™ RCS can be easily removed from the root canal if a retreatment is required.

**Clinical performance**

Permanent root canal obturation of mature teeth.

**Expected lifetime**

At least 5 years

**Mixing instructions**

- Aerate the powder by inverting the bottle of powder over twice and open the bottle.
- Using the spoon supplied in the box, take a leveled spoonful of powder and transfer it to the mixing pad.
- Detach a single-dose container of liquid. Twist the cap to open. Deposit 5 drops from the single-dose container on to the mixing pad without releasing the pressure.
- Gradually add the powder to the liquid and mix until obtaining a smooth paste (about 60 seconds). The mix is to be prepared immediately before use.
- Close the bottle of powder tightly and quickly after each use
- Immediately rinse and clean the instruments to remove any residual material.

**Placement into the root canal**

- Perform the root canal cleaning and shaping procedure using standard endodontic procedures.
- Select standardized gutta-percha master cone and check that it fits snugly at the working length.
- Flush with sterile saline solution to remove any remains of irrigant and dry without completely dehydrating the root canal
- Prepare BioRoot™ RCS according to the above mentioned mixing instructions.
- Apply a coating of BioRoot™ RCS onto the canal walls using either a paper point or the gutta-percha point.
- Complete the obturation by inserting the gutta-percha master cone previously coated with BioRoot™ RCS (single cone technique) or several coated gutta-percha points (lateral condensation technique).
- It is recommended to take X-rays to ensure that root canals have been thoroughly filled and that the working length does not exceed 2 mm from the apex.
- Use a hot instrument to sever the protruding gutta-percha.
- Perform the temporary restoration using any temporary filling material.
- If a post placement is required, it should be performed after complete hardening of BioRoot™ RCS, i.e. between 1 week and 1 month.

**Retreatment procedure**

Remove the existing sealer (BioRoot™ RCS or another sealer) using conventional removal techniques. Then follow instructions given in the sections "Mixing instructions" and "Placement into the root canal".

**Warnings and precautions for use**

It is at the discretion of the practitioner to choose the product according to the patient's clinical condition and current endodontic practices.

**Precautions for use**

Use the powder within 6 months after first opening.
Liquid in single-dose containers must be used immediately after opening.
Use the 5 liquid single-dose containers within 3 months after first opening of the foil.

**Warnings**

- Do not use the mixed paste on several patients to avoid cross contamination.
- The working length in the root canal should not exceed 2 mm from the radiological apex to avoid apical extrusion.
- The device may cause eye irritation and skin burn. Do not manipulate the device without appropriate personal protective equipment.
- Use a rubber dam to ensure that the cement is not in contact with the mucosa or swallowed.
- Use the spoon provided with BioRoot™ RCS. Otherwise, technical specifications are not guaranteed and may lead to the failure of the endodontic procedure.
- Do not use BioRoot™ RCS powder with another liquid than BioRoot™ RCS liquid. Otherwise, technical specifications are not guaranteed and may lead to the failure of the endodontic procedure.
- Use the recommended ratio of liquid/powder. Otherwise, technical specifications are not guaranteed and may lead to the failure of the endodontic procedure.

• Do not use a warm technique obturation. Otherwise, technical specifications are not guaranteed and may lead to the failure of the endodontic procedure.
• Using BioRoot™ RCS without gutta-percha points jeopardizes endodontic and retreatment procedures.

**Side effects**

Possible hypersensitivity to BioRoot™ RCS

**Disposal**

No specific recommendations

**Storage**

Do not store above 30 °C; protect from moisture.

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Assembly kit for 35 applications, containing:

- 35 single-dose containers of 0.25mL liquid divided into 7 foils
- 1 bottle of 15 g powder
- 1 measuring spoon

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